


OVERVIEW

During an Upper Endoscopy , an Urease (CLO) test may be taken. MEDILOGIK EMS™ enables the Endoscopist to record the location where the test was taken. Depending on the type of test kit used it is also possible to record the test result.

RECORDING A UREASE TEST

1. Proceed through the care pathway for the Upper Endoscopy procedure to the “Findings” screen.
2. On the anatomical diagram, mark the site where the test was taken. A numbered mark will appear.
3. To the right of the diagram, select the “Site” tab and observe the **[Urease Test]** button at the bottom of the screen.
4. Click the **[Urease Test]** button, to record that the test was taken, observe the dialogue window that opens.
5. Click the checkbox to confirm that the test has been taken.
6. If you know the results of the test you may record positive or negative, alternatively you can click **[Close]** to shut the dialogue window.
7. The remainder of the report is completed as normal.



RECORDING UREASE TEST RESULTS POST PROCEDURE

If your organisation uses Urease test kits that take some time to provide results, please follow these steps to record the findings. Users will be given the opportunity to reprint the procedure reports if required:

1. Log in to EMS.
2. Select the “Outstanding Urease Results” panel.
If this is not visible to you then you can enable it from the **HOME / PANELS** menu. Select “Outstanding Urease Results” and the panel will appear.
3. Locate the correct patient to be updated.
4. Click **[Record Result]** next to their name.
5. Confirm whether the result was positive or negative and click **[OK]**.
6. The system will then present the user with the option to reprint the updated procedure report.
7. Choose **[Yes]** to reprint and **[No]** to return to the previous screen.

