

Software Release Note – Version 1.10.4

Application	Version	Release Date	Pre-requisite
MEDILOGIK EMS	1.10.4	17 Nov 2016	<ul style="list-style-type: none"> MS Silverlight V5.1 .NET 4.6.1 on server

Important Notes

This version **requires** a minimum of .NET 4.6.1 to be installed on the server prior to upgrade. This version of .NET is backwards compatible with the previously required minimum of .NET 4.5.

New Feature

N/A

Improvement

- [EMS-1710] - Added additional SciStore export configuration settings to allow including hospital identifier, and to generate export files if CHI and/or Endoscopist code are missing. All options are in an off state by default therefore this will not affect existing SciStore exports.

(This provides an option for a private HealthCare customer to use this as a pdf report export)

Bug

- [EMS-1694] - If no scopes existed within the system for a given site & procedure type the "Add" scope button was incorrectly visible on the In Room step.
- [EMS-1695] - Crash when moving directly to Recovery step (very uncommon permission mix) or clicking Complete button (an unnecessary action) on the Report step of the first procedure of a multi-procedure.
- [EMS-1708] - **ABMU ONLY** - Fix to treat QRY/ADR interface as at ABUHB (with specific organisation code matches) so that NHS numbers are mapped correctly.

MEDILOGIK EMS Contact Details

Please address all support calls to:

Email: support@medilogik.co.uk

Tel: +44 (0)1473 351666