

## Software Release Note – Version 1.10

Application	Version	Release Date	Pre-requisite
MEDILOGIK EMS	1.10	9 <sup>th</sup> August 2016	MS Silverlight V5.1 .NET 4.6.1 on server

### Important Notes

This version **requires** a minimum of .NET 4.6.1 to be installed on the server prior to upgrade. This version of .NET is backwards compatible with the previously required minimum of .NET 4.5.

### New Feature

- [EMS-1486] - Provided instructional text for the user when they log in with no configured panels (typically on first login)
- [EMS-1543] - Redesigned the **InRoom screen**. The separate image tab for the live image feed has been removed and incorporated in the single InRoom screen.
- [EMS-1559] - Added support for integration with external scope tracking systems, including Fingerprint ScopeTrack and Getinge T-DOC
- [EMS-1592] - Added basic support for HL7 v2.4 ADT^A05 messages

### Improvement

- [EMS-1400] - Patient's previous history now automatically refreshes when performing navigation and save actions without having to reload the patient record
- [EMS-1518] - EMS Follow-ups panel performance improved and site filter added
- [EMS-1526] - The default date range selection for follow-ups / cancellations is now weeks, not days
- [EMS-1527] / [EMS-1570] - Various National Endoscopy Database (NED) configuration improvement to assist in pilot testing with first customer.
- [EMS-1542] - Improved the layout of the Session List - Outcomes report to avoid unnecessary wrapping

### Bug

- [EMS-1447] - Scheduling: It should not have been possible to create a case mix with both EMS and Non-EMS procedure types
- [EMS-1477] / [EMS-1564] / [EMS-1565] - Letter editor: Fixed some merge fields to include/exclude the year in the date correctly
- [EMS-1488] - Fixed various styling and alignment issues
- [EMS-1498] - When editing a patient manually the CHI number became 10 space characters if the number was previously removed. This field is now correctly blanked when emptied
- [EMS-1541] - Scheduling: Previous templates no longer show for deleted Endoscopists
- [EMS-1558] - The visibility of "Edit Patient" and "Merge This Patient Into Another Patient" in the Patient menu should respect the global configuration setting "Can Manually Add Patients", but didn't (the options performed no action when selected, but they should not have been visible)
- [EMS-1568] - The Summary Panel sometimes incorrectly displayed more than one site even when only a single site was selected

## Software Release Note – Version 1.10

- [EMS-1574] - The Summary Panel sometimes did not include the full count of cancellations recorded on the day
- [EMS-1577] - The Completion Confirmed by field on the Summary screen had incorrect mandatory asterisk marking
- [EMS-1581] / [EMS-1582] - The incoming referrals and deaths and readmissions panels details windows had incorrect mandatory validation asterisks and some alignment issues
- [EMS-1586] - Bronchoscopy Post-Pathology Diagnosis of Foreign body had incorrect qualifiers
- [EMS-1588] - Selecting only a Trainer in room would cause export to SCI Store to fail

### MEDILOGIK EMS Contact Details

Please address all support calls to:

Email: [support@medilogik.co.uk](mailto:support@medilogik.co.uk)

Tel: +44 (0)1473 351666