

Software Release Note – Version 1.9

Application	Version	Release Date	Pre-Requisite
MEDILOGIK EMS	1.9	04-May-2016	MS Silverlight V5.1 .NET 4.6.1 on server

Important Notes

This version **requires** a minimum of .NET 4.6.1 to be installed on the server prior to upgrade. This version of .NET is backwards compatible with the previously required minimum of .NET 4.5.

System Administration/Configuration

- [EMS-1375] – Permission to Add/Modify Staff restricted to System & Site Administrator roles
- [EMS-1414] - Config setting to default copies of Reg and Ref GP report copies to zero if enabled to despatch via EDT
- [EMS-1428] - Config setting to allow entire Patient History to be automatically visible when loading Assessment step
- [EMS-1430] - Removed buttons on Referral step for "New GP" and "New Consultant". This functionality should be used from the scheduling application
- [EMS-1397] - Added ability to provide custom messages to inform users of requirement to contact their own helpdesk for failed logins

Interfacing/Interoperability Enhancements

- [EMS-1339] - Implemented export to Fairwarning system
- [EMS-1458] - Implemented MEDITECH patient demographics interface
- [EMS-1399] - Created "Docman EDT Issues" dashboard management panel
- [EMS-1410] - Docman EDT - Modify ReportType to include actual procedure name
- [EMS-1478] - Docman EDT - Vary Department value by procedure type

General Performance Improvements

- [EMS-1347] - OGD Performance Report
- [EMS-1349] - Referral By Appointment Type Report
- [EMS-1357] - Colonoscopy Performance Report
- [EMS-1364] - EMS follow-ups panel
- [EMS-1365] - Waiting List Panel
- [EMS-1366] - Awaiting Vetting Panel
- [EMS-1368] - Cancellations Panel

New Feature

- [EMS-1336] – Mandatory fields indicated with a red asterisk (*) throughout the application

Software Release Note – Version 1.9

Scheduling Improvements

- [EMS-1383] - Scheduling grid includes a separate site row above room name
- [EMS-1398] - Scheduling grid has the room and session period fixed in place above the grid, allowing to scroll without losing sight of the room etc
- [EMS-1206] - Users can select all sites in cancellations follow-up, not their allocated sites
- [EMS-1382] - Summary Panel content enhanced and performance improved
- [EMS-1420] - Session List renamed to “**Session List– Scheduled**”
- [EMS-1420] - Booked points has been added to the “Session List– Scheduled”
- [EMS-1423] - Single or multiple list options when viewing “Session List– Scheduled”
- [EMS-1418] - New report “**Session List - Outcomes**” to enable Endoscopist to manually record OPCS4 and Outcome Codes
- [EMS-1279] - New report “Procedure Numbers” available by Unit and date detailing Endoscopist and Procedure Types

Clinical Reporting Improvements (in pathway step order)

- [EMS-1411] - **Patient History** - added ability to go direct to view images (for use at MDT)
- [EMS-1409] - **Daily List Panel** - points allocated to procedures now shown
- [EMS-1429] - When loading an episode from the Daily List panel if Referral Reasons haven't been recorded the episode will be loaded directly at the Assessment step enabling the Endoscopist to view High Risk issues and current medication
- [EMS-1387] - **Findings** - Ability to record Urease test on ERCP and EUS procedures
- [EMS-1431] - **Summary** - Intended & Actual Extent now displayed in logical anatomical order
- [EMS-1415] - **Summary** - Hypoxia added to Procedure Complications for all procedures
- [EMS-1297] - **Diagnosis** - Extrinsic Compression added to OGD diagnosis
- [EMS-1388] - **Diagnosis** - single “Normal” button added for OGD diagnosis
- [EMS-1412] - **Post Procedure Diagnosis** – Added Coeliac disease (OGD, EUS and ERCP)

Procedure Report

- [EMS-1371] - Changed “GP” to read “Registered GP” taken from Referral Details
- [EMS-1422] - Report print option includes Trainee if present
- [EMS-1342] - Report includes a list of report recipients below Endoscopists name
- [EMS-1378] - The Images Report in a combined procedure include “No Images Were Captured” when one of the procedures has no images saved
- [EMS-1392] - The Histology Request in a combined procedure states no biopsies have been taken when only one of the procedures has pathology recorded

Software Release Note – Version 1.9

Bug Fixes

- [EMS-1353] - Fixed minor spelling mistake on Summary panel
- [EMS-1358] - When scheduling an appointment, the Endoscopist dropdown did not update on calendar month changes
- [EMS-1359] - Fixed prompt to save on navigate occurring twice
- [EMS-1363] - Improved memory usage of Silverlight client by fixing various memory leaks
- [EMS-1424] - OGD Performance Report - Mean Analgesia drilldown caused server error
- [EMS-1444] - Unscheduled/Emergency procedures are being incorrectly shown in the Scheduling Activity Report
- [EMS-1465] - Cancellation and Rejected Referrals report final drilldown did not show any data when clicking specific cancellation reasons
- [EMS-1473] - Hide site filter on Summary panel when user has only a single site in their profile
- [EMS-1348] - Addition of Image Taken for D2 Completion on OGD Perf Report

MEDILOGIK EMS Contact Details
Please address all support calls to support@medilogik.co.uk Support Line 01473 351666